# Vær oppmerksom på hudflenger

Forstå og håndter hudflenger

Lett

å

overse

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**Referanser:** 1. Malone, L.M., et al. The Epidemiology of Skin Tears in the Institutionalized Elderly. *Journal of the American Geriatrics Society*, 1991. 2. Konya C. et al. Skin injuries caused by medical adhesive tape in older people and associated factors. *Journal of Clinical Nursing*, 2010. 3. Carville K. et al. STAR: a consensus for skin tear classification. *Primary Intention: the Australian Journal of Wound Care*, 2007. 4. LeBlanc K. et al. Skin tears: The forgotten wound. *Nursing Management*, 2014. 5. Bryant R. A. Types of Skin Damage and Differential Diagnosis In Bryant, R.A., Nix D.P., *Acute and Chronic Wounds: Current Management Concepts*, edition 5, St. Louis, Elsevier, 2016. 6. LeBlanc K. et al. Best practice recommendations for the prevention and management of skin tears in aged skin. *Wounds International*, 2018. Nedlastbar fra [www.woundsinternational.com](http://www.woundsinternational.com). 7. White MW. et al. Skin tears in frail elders: A practical approach to prevention. *Geriatric Nursing*, 1994. 8. Resnick B. Wound care for the elderly. *Geriatric Nursing*, 1993. 9. LeBlanc K. et al. International skin tear advisory panel: A tool kit to aid in the prevention, assessment, and treatment of skin tears using a simplified classification system. *Advances in Skin and Wound Care*, 2013. 10. Woo K. & LeBlanc K. A randomised controlled pragmatic study to evaluate the use of silicone dressings for the treatment of skin tears. Poster presentation at Wound Con Summer [virtual] conference, 2020. 11. Patton P. et al. An open, prospective, randomized pilot investigation evaluating pain with the use of a soft silicone wound contact layer vs bridal veil and staples on split thickness skin grafts as a primary dressing. *Journal of Burn Care and Research*, 2013. 12. David F. et al. A randomised, controlled, non-inferiority trial comparing the performance of a soft silicone-coated wound contact layer (Mepitel One) with a lipidocolloid wound contact layer (UrigoTul) in the treatment of acute wounds. *International Wound Journal*, 2017. 13. Collin O. Use of Mepitel One dressing following hand surgery: a case study series. Poster presentation at Wounds UK Conference, United Kingdom, 2009. 14. Edwards J. et al. Hand burn management: minimizing pain and trauma at dressing change. *British Journal of Nursing*, 2013. 15. White R. A multinational survey of the assessment of pain when removing dressings. *Wounds UK* 2008. 16. Nelson D. Better outcomes for skin tears with new 5 layer bordered foam dressings. Poster presentation at 50th Annual Conference of the Wound Ostomy and Continence Nurses' Society, Philadelphia, Pennsylvania, United States of America, 2018. Quality Improvement Project. 17. Rook S, Davies P, Frenthoff E, Würfel T. Mepilex® Border Flex — results of an observational study in German specialist wound care centres. *Wounds International*. 2019;10(1):40-3. 18. Data on file 19. Meaume S. et al. A study to compare a new self-adherent soft silicone dressing with a self-adherent polymer dressing in stage II pressure ulcers. *Ostomy Wound Management*, 2003. 20. Woo K.Y. et al. A randomized crossover investigation of pain at dressing change comparing 2 foam dressings. *Advances in Skin and Wound Care*, 2009. 21. Stephen-Haynes J. et al. The assessment and management of skin tears in care homes. *Br J Nurs*, 2011. 22. Gray D. et al. Pilot RCT of two dressing regimens for the management of skin tears. *Wounds UK*, 2011.